STEREOTACTIC BREAST BIOPSY DIGITAL REPEAT ANALYSIS SUMMARY

Site: _			
Digita	Stereotactic Breast Biopsy System:		
Date	From: To:		<u></u>
	Cause	Number of Exposures	Percentage of Repeats
1.	Positioning		
2.	Patient Motion		
3.	Detector Underexposure (excessively noisy images)		
4.	Improper Detector Exposure (saturation)		
5.	Incorrect Patient ID		
6.	X-Ray Equipment Failure		
7.	Software Failure		
8.	Blank Image		
9.	Other		
Repe Rat			
	al Repeated SBB Films 1-9 al SBB Films Used		
100	Repeat Rate		
Perfor	med by	Date:	
	for Quality Improvement		
Physician Review		Date:	
Medical Physicist Review		Date:	

Figure 10. Digital stereotactic breast biopsy repeat analysis summary